

MRI Simulator

Objective:

This tool will be used to simulate roughly the steps of operating an MRI scanner. This includes planning of the exam, technical considerations, and the additional component of grading and seeing a final image.

I am not sure which would be easiest or most cost effective to host this simulation. Would it be easier to host:

- As a standalone activity?
- Wordpress plugin?

Step 1: Planning	Step 2: Technical Considerations	Step 3: Grading
<p>This includes:</p> <p>There should be 3 viewers with different sets of anatomical images within them. The boxes will be labeled anatomically as well and should communicate with each other.</p> <p>Defining two areas:</p> <ul style="list-style-type: none">• 1 area that needs to be covered. This is an area of the image that requires users to plan a coverage area through it.• 1 area that defines areas that need to be accounted for. This means that users need to make sure that areas outside the first box are covered as well, but do not require users to plan through it. <p>We also need to define a correct angle. This will relate to the other 3 viewers.</p>	<p>This includes:</p> <p>A list of parameters that can be changed by the user. The changes will:</p> <ul style="list-style-type: none">• Change the area of coverage in the viewers.• Produce calculations that will result in a grade.	<p>A list of parameters and a resulting score that correlates with it that will add up to a final score. This score should be able to be exported. Along with a grade, a stock photo will be displayed demonstrating a final image. Image filtering should manipulate the image based on the calculation results.</p> <p>Ex.</p> <ul style="list-style-type: none">• <i>Parameter SNR will apply a noise filter on the image making it look grainy.</i>• <i>Parameter CNR will apply a blur filter on the image to make it look blurry.</i>• <i>Parameter resolution will apply a sharp filter to make things look more detailed.</i>• <i>Parameter FOV will apply a zoom filter to make thinks zoom in or out.</i>

Terminology

Activity: This is an entire lesson for a user, it contains individual simulations

Simulation: An independently graded activity situation a user will be challenged with

Learn Mode: This is a progression where users are guided to the correct answer.

Test Mode: This is a progression where users are tested on their performance.

Template: This is a format where the viewer and parameters will be laid out.

Viewer: This is a image box that contains multiple images. A user will work at covering the correct area on these images.

Parameters: A display of values that a user must manipulate.

High Level Parameters: The user will change parameters which will alter the high-level parameters. The user can not directly change the high-level parameters.

RF Coils: This is a selection by the user that will determine the results of their performance.

Field Strength: This is a selection defined by the admin on activity set up which will alter the result of the high-level parameters.

Question Exercise: This is an exercise used in place of a simulation. It will ask a question, teach a concept, or present information.

Parameter Exercise: This is an exercise that will test the user's ability to manipulate parameters and perform an MRI.

Anatomy Laterality: This is meant to describe the sides of our patient. (Right (R), Left (L), Anterior (A), Posterior (P), Foot (F), Head (H))

Image Orientation: This describes how we slice our patient. (Sag = separate patient into right and left side) (Cor = Separate patient into anterior and posterior side) (Ax = separates patient into a head and foot side)

Question Exercise: This is a simulation type that will ask the user a question, get an answer, and possibly show a solution as to why the answer is right or wrong. We can also determine if we want this exercise seen in the simulation list area.

Simulation list area: This is the area where a user can see all the simulations they still need to complete.

Image Contrast: Based on values the user selects, the image will look different. These are always the same based on specific values.

Pulse sequence: This is a way we can perform a simulation. It has specific parameters, limits, and values associated with it.

R-R: When using the cardiac exercise, we based all our scan time values off this period of time. It represents the cardiac cycle.

Grade Card: This is a pop up of what the user scored after each simulation.

Activity Results: This is a page seen after an activity which gives an overall score.

Body Part: This section is used for tagging activities and eventually requiring specific body parts to be completed for a user.

Lesson Area: This is a section on the frontend where education to the user can be presented. This area can present lessons or other content.

Vendor to Vendor: Parameters can have different layouts (vendor = template). When a vendor to vendor activity is selected, to parameter selections can be view (for each template selected) stacked on top of each other.

Box 1: This is a box representing an area that needs to be covered by the user.

Box 1 angle: This is the angle of the box that needs to match within a couple degrees of the correct answer.





Box 2: This is a box representing an area that the user does not need to include with parameter FOV, but still account for.

Image Filtering: This is a stage of the process where high-level parameters control the appearance of the final image.

Global Setting

Each activity will include a collection of separate simulations. Example:

- Activity 1 (With Global Setting)
 - Simulation 1 > show image
 - Simulation 2 > show image
 - Simulation 3 > show image
 - Question 1
 - ...
 - Grade with score average

2 Modes of Operation	Templates (Optional)
<ul style="list-style-type: none"> • Learning mode: This will show descriptions and graphs demonstrating objectives. This will guide users through the planning of the exam and provide lessons. • Testing mode: This will be an open simulation that will present results for grading. <p>In both modes, the option to provide a simulation or question slide within an activity should be an option.</p> <div> <div> Parameters TR TE TI ETL etc. (Expanded on later) </div> <div> Questions Title Image/video/shortcode Question/Statement Answer Answer Answer Answer Solution </div> </div>	<p>This provides different layouts that I can choose to build with. In other words, the layout of the parameters will be different based on one of 4 templates I choose. All the parameters will have the same purpose, but will be displayed differently.</p> <ul style="list-style-type: none"> • RITE Advantage • Siemens • GE • Toshiba • Phillips <div>   </div> <div>   </div>

Backend Global Setting for Activity

- Activity Name:
- Instructions: (To explain the purpose of the activity)
 - Text area where instructions can be provided to user
- Mode: (To define how the simulations will be handled and what the goal for the user is)
 - **Learn**: Run stimulation in a learning setting defined below.
 - **Test**: Test Users
- Template: (Optional parameter display)
 - **Siemens**: Display parameters with a specific template
 - **GE**: Display parameters with a specific template
 - **Phillips**: Display parameters with a specific template
 - **Toshiba**: Display parameters with a specific template
 - **RITE Advantage**: Display parameters with a specific template
- RF Coil: (Used to test user on correct setup and for high-level parameter calculation)
 - Body
 - Knee
 - Brain and Neck
 - Brain T-R
 - Surface Coil
 - Foot
 - Breast
 - Shoulder
- Field Strength: (Used for high-level calculation)
 - 1.5T
 - 3T
- Body Part: (Used for advanced activity tracking)
 - Brain
 - IAC
 - Pituitary
 - MRA head
 - MRV head
 - Orbit
 - Neck
 - MRA neck
 - Cervical spine
 - Thoracic spine
 - Lumbar spine
 - Sacral spine
 - Chest
 - Breast
 - Cardiac
 - Pelvis Boney
 - Pelvis Soft Tissue
 - Hip
 - Leg Bone
 - Knee
 - Ankle/Foot
 - Shoulder
 - Elbow

- Wrist
- Hand

Backend Set-up

Activity Title: _____		Field Strength: <input type="radio"/> 1.5T <input type="radio"/> 3.0T	
Title: _____			
Instructions: _____ _____		Acceptable RF Coils:	
Modes:		<input type="checkbox"/> Whole Body Receiver <input type="checkbox"/> Knee (Transmit-Receive) <input type="checkbox"/> Brain/Neck (Receive only) <input type="checkbox"/> Brain (Transmit-Receive) <input type="checkbox"/> Surface Coil (Receive only) <input type="checkbox"/> Foot/Ankle Coil (Receive Only) <input type="checkbox"/> Breast (Receive only) <input type="checkbox"/> Shoulder (Receive only)	
<input type="radio"/> Learn	<input type="radio"/> Test	Vendor	
		<input type="text"/> RITE Advantage GE Siemens Phillips Toshiba	
Upload Lesson or Enter Shortcode	Upload Lesson or Enter Shortcode	or	
		<input type="text"/> to <input type="text"/> RITE Advantage GE Siemens Phillips Toshiba	
+ <input type="text"/> Question Exercise Parameter Exercise		Body Part:	
		<div style="display: flex; justify-content: space-between;"> <div> <input type="radio"/> Brain <input type="radio"/> IAC <input type="radio"/> Pituitary <input type="radio"/> MRA head <input type="radio"/> MRV head <input type="radio"/> Orbit <input type="radio"/> Neck <input type="radio"/> MRA neck <input type="radio"/> Cervical spine <input type="radio"/> Thoracic spine <input type="radio"/> Lumbar spine <input type="radio"/> Sacral spine <input type="radio"/> Chest <input type="radio"/> Breast </div> <div> <input type="radio"/> Cardiac <input type="radio"/> Pelvis Boney <input type="radio"/> Pelvis Soft Tissue <input type="radio"/> Hip <input type="radio"/> Leg Bone <input type="radio"/> Knee <input type="radio"/> Ankle/Foot <input type="radio"/> Shoulder <input type="radio"/> Elbow <input type="radio"/> Wrist <input type="radio"/> Hand </div> </div>	

Frontend View


Title: _____	Body Part:	Field Strength:
Instructions:	Select RF Coil:	
<input style="width: 100%;" type="text"/>	<input type="checkbox"/> Whole Body Receiver <input type="checkbox"/> Knee (Transmit-Receive) <input type="checkbox"/> Brain/Neck (Receive only) <input type="checkbox"/> Brain (Transmit-Receive) <input type="checkbox"/> Surface Coil (Receive only) <input type="checkbox"/> Foot/Ankle Coil (Receive Only) <input type="checkbox"/> Breast (Receive only) <input type="checkbox"/> Shoulder (Receive only)	

Step 1

When setting up the simulation, a specify way to test a user will be displayed. There are two paths.

1. A list of parameters that can be manipulated and graded.
2. A popup slide that can be shown on the frontend or hidden that present messages or questions.


On the backend, the admin should be able to select what they want to add.



If Question Exercise is selection see part 2. No image planning is required for this type.

If Parameter Exercise is selected, continue to image planning.

Simulation Builder: Image Planning	
Backend Builder	Frontend User view
<ul style="list-style-type: none"> Area to upload multiple images into one of 3 separate viewers Two boxes <ul style="list-style-type: none"> Box 1 = Area that needs to be covered with the following parameters (<i>See next section</i>) <ul style="list-style-type: none"> FOV Slice thickness Gap Box 2 = Area that needs to be accounted for with the following parameters. <ul style="list-style-type: none"> FOV RectFOV Anti-Aliasing Box 1 angle <ul style="list-style-type: none"> I would like to set a correct angle on the image with a degree of acceptance. This will give a range of correct answers. The angle will correlate with a plane selected Axial (parallel to R-L and A-P, but perpendicular to H-F) Coronal (parallel to R-L and H-F, but perpendicular to A-P) Sagittal (parallel to H-F and A-P, but perpendicular to R-L) Other = don't grade <p>Rough view for Backend Builder</p> <p>Each of these viewers can have 3 to 4 images of the same body part. I should be able to load these images.</p>	<p>The frontend user will be given a box to manipulate. They can move the box anywhere on the image and angle any way they please.</p> <p>Each viewer allows for images to be paged through They should be able to click through these images to make sure they are positioned correctly.</p> <div> </div> <p>Learn mode= If the user angles different than how it was built, the box will turn red. This will not let user continue until they have the correct angle</p> <p>Rough View Learn Mode</p> <div> </div> <p>Test mode = No sign of correct or not will be given, it will be graded as correct or incorrect. This view won't turn the box red if the position or angle is off.</p>




The **yellow box** represents “Box 1”. This tests the user on the angle and area of coverage. This is the only box the user sees.

The **green box** represents “Box 2”. This defines the rest of the anatomy that needs accounted for in the image.



Image Labeling and Function


Backend Design Idea



“Simulation Title” will be displayed in the “simulation view”. The “image label” is the gallery descriptor (*In the image below, it is “Sagittal T1”*).

The “anatomy laterality” is defined by the image orientation.







Then we should be able to set the boxes. “Box 1” sets the area that needs to be covered. We can refer to this as the **Field of View**.

The yellow line on the image on defines the slice direction.


Sag slice =




Cor Slice =



Axial Slice=





Then box 2 can be placed.




The slice direction will remain as defined (sagittal, coronal, axial) unless the angle is ± 45 degrees. Then it will default to the next slice direction.

In the image below, it is planned as a sagittal. If we move beyond ± 45 degrees, it becomes a coronal. (see image planes above for slice directions).



The slice direction is defined by *slice thickness* and *gap*. (see next section). As we increase the slice thickness and gap, we increase the coverage in the slice direction and visa versa.

Additional Labelling



Any move of the field of view will reflect across all viewers. Example: If the field of view is moved to the "R", it will move to the "R" on the axial and coronal. If it moves towards the "A" it will move on the axial and sagittal viewers. Etc..
(See Video)

Step 2: Technical Considerations

A list of parameters will be found on the template that has been chosen. There will be two types of simulations.

1. Question exercise
2. Parameter exercise

1. If question is selected, information would be required to display appropriately on the frontend.

Title Message: _____

Upload MP4 or Picture
or
Imbed iframe with shortcode: _____

Question: _____
Answer: _____
Answer: _____
Answer: _____
Solution: _____
☐ Hide of Frontend
☐ Grade at end ☐ Retry Until Pass

Do not show in simulation list:

Simulation 1
Simulation 2
Question 1
Question 2
...

Let user answer then move on.

Make user try until they get it correct.

It would be nice to have the ability to also post shortcode for other plugins (such as H5P) that can be presented in an iframe.

2. If Parameter Exercise is selected, a list of parameters should be displayed. The option to input default fields would be needed.

Hide

Parameter: Default Value

High-Level Parameters

Limits

Hide and Don't Grade

☐ Slice Thickness ____mm
☐ Slice Gap ____mm
☐ FOV ____mm
☐ PE ____
☐ FE ____
☐ # Slices ____
☐ Phase Direction ____
☐ Alias On | Off
☐ RectFOV ____%
☐ NSA ____
☐ NA ____
☐ Saturation FS | WS
☐ Restore On | Off
☐ Half Fourier ____%
☐ Parallel Imaging On | Off
☐ Confirm Frequency Yes | No (use a checkbox)
☐ Cardiac Gating Retrospective | Prospective
☐ Cardiac Segment ____
☐ Cardiac Phase ____
☐ Trigger Window ____ms
☐ Trigger Delay ____ms
☐ Acquisition Window ____ms to ____ms (answer in between)
☐ Measurement ____
☐ RBW ____Hz/px
☐ ETL ____
☐ Gradient Mode On | Off
☐ RF Mode On | Off
☐ VENC ____cm/s
☐ B Value ____
☐ R-R ____ms
☐ Breathhold On | Off

Only used for grading purposes

Image Contrast

T1

T2

PD

IR: STIR

IR: FLAIR

Phase Contrast

DWI

TOF

SSFP

load

load

load

load

load

load

load

load

Correct image

Incorrect image

Pulse Sequence

SE

GRE

DWI

Phase Contrast


TOF

SSFP

Only used for grading and to auto hide parameters

+

Question Exercise
Parameter Exercise

Here is a list of the parameter and functionality. If the “learn mode” is selected, an  can be found next to the parameter and when clicked, the “description” will be displayed.

Parameters						
Parameter	Abbreviation	Min	Max	Purpose	Description	Field display
TR	TR	$((2 * 1 + [ES]) + ([FE] * (1 / ([RBW] / 2))) * [Num\ slices] * [ETL]) + ([fs] * ([TR] / [ETL] / 2)) / [num\ acquisition])$	10,000	Use of grading	Repetition time defines the period from one excitation pulse to the next. It refers to the cycle of repeating patterns of data collection. TR controls T1 contrast in our image. As we increase TR, we increase inherent signal to noise in our image, decrease T1 contrast, and increase scan time.	Number
TE	TE	$(2 * 1 + [ES]) + ([FE] * (1 / ([RBW] / 2)))$	1000	Use of grading	Echo time defines the period from excitation to when our echo is collected. Most pulse sequences demonstrate our TE within our TR. This means that our TR will always be longer than our TE in most situations. TE controls T2 contrast. As we increase TE, we increase T2 contrast, decrease inherent signal to noise, and increase susceptibility in our image.	Number
TI	TI	-	-	Use of grading	Inverse time is used to suppress tissues based on their T1 recovery. This means that inverse sequences influence T1 contrast. By utilizing a 180-degree RF pulse prior to our TR, we can invert protons within a volume of our patient. If we wait a period of time that correlates with the T1 of a tissue we wish to suppress prior to excitation, that tissue will be removed from our image. Shorter inverse times will correlate with fat suppression. We call these short T1 inverse recovery sequences or STIR. Longer inverse times may correlate with cerebral spinal fluid, and these are called fluid attenuation inverse recovery sequences or FLAIR.	Drop Down If 1.5T show 160 2500 If 3T show 220 2500
Flip Angle	FA	If gradient echo=10 If spin echo =120	If gradient echo = 90 If spin echo = 180	Use of grading	A flip angle defines the amount of energy we deliver to our patient to tilt hydrogen away from longitudinal magnetization. A spin echo sequence always starts with a 90-degree RF pulse. It then is followed by one or more 180 degree refocusing pulses. Most of the time, we can change this number to lower patient heating. We should not decrease it below 120 degrees however Gradient echo sequences allow us to use a flip angle of 90 degrees or less. As we decrease are flip angle, we decrease T1 contrast in our image. We also decrease patient heating.	Number

Slice Thickness	ST	0.1mm	50mm	Increase the coverage in the “slice direction”	Slice thickness controls how much tissue we will average into a 2-dimensional slice. The thicker our slice thickness, the more tissue we are averaging into our image. This can lead to image blurring called partial average volume artifact. Increasing our slice thickness can also promote higher SNR and increase coverage. As we increase coverage, we may be able to use less slices which leads to shorter scan times.	Number
Slice Gap	SG	0	5	Increase the coverage in the “slice direction”	Slice gap is a parameter that controls the spacing in between slices. The main benefit of this is to increase coverage so that we need less slices. This leads to shorter scan times. The issue with utilizing a slice gap is that pathology may be missed.	Number
Field of View	FOV	5	500	Increase the size of the “field of view” in viewer	The field of view controls the coverage included in our image. It acts as a frame around our image. Based on other parameters, we will have a fixed number of pixels in our image. As we decrease our field of view, we squeeze those pixels together making them smaller. If our pixels get smaller, our spatial resolution increases. This means that a decrease in field of view leads to an increase in spatial resolution. The downside of this is that our signal to noise ratio will decrease. This is because there is an inverse relationship between spatial resolution and SNR.	Number
Phase Encoding	PE	–	–	Use of grading	Phase encoding defines the number of lines in our data space or K space. As we increase our phase encoding, we increased the number of lines that need to be filled with echoes. This means we need to produce more echoes which leads to an increase in scan time. As we increase the number of lines however, we increase our spatial resolution. This parameter combined with our frequency encoding direction make up our image matrix or a grid of pixels which represent our patient. As we increase the number of phase encoding lines in our image, they are squeezed together due to a fixed field of view, and this leads to increased spatial resolution. There is an inverse relationship between spatial resolution and SNR which means that as we increase our phase encoding, we decrease our SNR.	Drop Down 128 192 224 256 320 384 448 512
Frequency Encoding	FE	–	–	Use of grading	Frequency encoding is a parameter that defines how many data points we will use to represent patient frequencies in the frequency encoding direction. As we increase our frequency encoding, we increase the number of data points used to replicate our patient signal. This can increase our echo spacing however it does not always increase our scan time. As we increase our frequency encoding, we increase the number of lines in our frequency encoding direction squeezing them together due to a fixed field of view. This means that we will increase our spatial resolution which leads to a decrease in SNR due too in inverse relationship between the two.	Drop Down 128 192 224 256 320 384 448 512
Number of Slices	#Slice		80	Increase the coverage in the “slice direction”	The number of slices define how many pictures we need to cover in area of our patient. This is limited to how long our repetition time is and how many echoes we are collecting per repetition time.	Number
Phase Direction	Phase Dir	–	–	Determine if “Box 2” is covered properly. The area represented by the green box needs to be accounted for.	Phase encoding direction defines which orientation our phase encoding lines will be filled. This is important due to the many artifacts seen in this direction. The big artifacts to worry about in our phase direction is aliasing or wrap artifact and motion.	Drop Down A-P (Show if Sagittal or Axial) H-F (Show if Sagittal or Coronal) R-L (Show if Coronal or Axial)
Anti-Aliasing	Alias	–	–		Any tissue outside of our field of view in the phase encoding direction can lead to aliasing artifact. This parameter will compensate for this accounting for tissue outside our field of view in the phase encoding direction.	
RectFOV	RectFOV	–	–	This will increase or shorten the field of view in the “phase direction” by the percentage of our FOV defined in the drop down	Rectangular field of view is a parameter that will add or subtract information collected in the phase encoding direction. Is important to remember that we are not changing the number of lines in our k space and therefore this parameter will not influence our spatial resolution. It will however alter how much signal we collect from our patient and therefore as we decrease our rectangular field of view, we will decrease our SNR. Also, if we decrease our rectangular field of view, we no longer need to collect as many lines of k space to build an image and therefore can reduce our scan time.	Drop Down 200% 175% 150% 125% 100% 93% 85% 75% 68%
Number of Signal Averages	NSA	1	10	Use of grading	There are many situations where we have a signal deprived image. The number of signal averages can be increased so that we can increase our signal to noise ratio. What happens when we do this is that we collect multiple sets of information for the same slice depending on our signal averages we input. We can then average these data sets together to	Number

					produce in average data set. Since noise is random, we can decrease the noise in our image by this factor $\sqrt{2}$.	
Number of Acquisitions	NA	1	10	Use of grading	The number of acquisitions is a parameter that allows us to separate the echoes collected for a given TR into multiple measurements. The result is that we can use a shorter repetition time or utilize multiple breath holds when needed.	Number
Saturation	S	–	–	Use of grading	It is important to suppressed tissue in our image too optimize pathological processes. We can suppress fat, water, or silicone based on the unique frequencies They have. This is very sensitive to the homogeny or uniformity of our magnetic field at isocenter. Because we're not dependent on tissue relaxation to suppress tissues in our image, we can use this technique pre and post contrast.	Drop Down Fat Sat. or Water Sat. = 1 Off = 0
Restore	R	–	–	Use of grading		Check box
Half Fourier	Half Four	–	–	Use of grading	Half Fourier is a technique that relies on the fact that our k space is a mirror image of itself. This means that the top half of our k space is a mirror image of the bottom half, and the right half is a mirror of the left half. Because of this, if we fill a minimum of 53% of our k space with echoes and fill the rest with the zeros, we can produce a diagnostic image. By doing this, we can reduce the number of echoes we need to collect and therefore decrease our scan time. Since we are collecting less echoes from our patient, our SNR will decrease when implementing this technique.	Drop Down 100% 87% 75% 63%
Parallel imaging	PI	–	–		Parallel imaging is a technique that can be used when we have two or more coil elements in the phase encoding direction. By using each coil element separately to collect data, we can produce multiple aliased images in a shorter scan time. By combining these datasets together, we can produce one full image in a shorter duration of time. This technique will shorten our scan time but also decrease our SNR.	Number On = 2 Off = 1
Confirm Frequency	Con Freq	–	–		Due to electron sharing throughout the electron cloud of fat molecules, it produces a slightly different magnetism than water and therefore precesses different than water. This produces chemical shift between fat and water. Confirming our frequency allows us to isolate a frequency we want to focus on to either suppress a tissue or enhance a tissue.	Check Box
Cardiac Gating	Car Gating	–	–	Use of grading	There are two ways we can perform cardiac gating. Prospective gating will image during diastole. Retrospective gating will image throughout the entire cardiac cycle and produce images at different stages along it. This leads to a dataset which looks like a movie of the heart pumping.	Check Box Retrospective (CINE) Prospective
Segments	Seg	1	20	Use of grading	When performing retrospective cardiac gating, we have the option to select how many points along the cardiac cycle we want to produce an image in. In other words, how many echoes we want to collect throughout the cardiac cycle. Increasing this number will allow for more images of the cardiac cycle however it will lead to an increase in temporal resolution which can cause artifacts in our image.	Number
Calculated Phases	Car Cal Phas	1	40	Use of grading	This parameter will allow us to reconstruct additional phases of the cardiac cycle based on phases we've already collected.	Number
Trigger window	Trig W	10ms	500ms	Use of grading	When performing prospective gating, we need to identify an acquisition window which is optimal for imaging. We use our R-R peaks to trigger imaging. In order to isolate our acquisition window, we must add a delay before our first R peak. This is our trigger window.	Number
Trigger Delay	Trig D	10ms	500ms	Use of grading	In order to appropriately identify or acquisition window, we utilize our R-R peak. In order to capture the point of our cardiac cycle where our heart is not moving for the longest duration of time, we must add a delay after our first peak. We call this the trigger delay.	Number
Acquisition window	Acq Win	$[R-R] - ([Trig D] + [Trig W])$		Use of grading	The period of time we collect data when using perspective gating is called our acquisition window. We identified this time using our trigger window and trigger delay. This period of time will always be less then the time it takes from one peak to the next. This causes an issue when trying to obtain T2 weighted images. When imaging T2 weighted images, we sometimes need to utilize more than one R-R peak.	Number
Measurements	Measure	1	100	Use of grading	The number of measurements defines how many times we will capture the same image per a given acquisition. This is utilized when performing a dynamic imaging. Temporal resolution becomes very important when we choose to use multiple measurements.	Number
Receiving Bandwidth	RBW	10hz/px	1700hz/px	Use of grading	Receiving bandwidth defines the rate that we sample our patient frequency at TE. Receiving bandwidth is tied closely with TE and our frequency encoding. As we increase our receiving bandwidth, we decrease our sample time which leads to possibly having the option to use shorter TR values. It also allows us to replicate our patient frequency more accurately and therefore decrease metal artifacts and blurring in our image. By decreasing our receiving bandwidth, we increase our sample time which increases our SNR and decreases our SAR. By doing this however we will decrease the contrast to noise ratio in our image.	Number
Echo Spacing	ES	$((1/\text{bandwidth}/\text{by acceleration factor (Parallel imaging)}) * [GM] * [RFM]) * 100$		Use of grading	Echo spacing refers to the period of time between echoes collected in a given repetition time. When multiple echoes are collected during a TR, there is a delay in between them. This is called echo spacing.	Only Show
Echo Train Length	ETL	1	[Pe]	Use of grading	Echo train length is a parameter utilized when performing spin echo sequences. It defines how many echoes we will collect within a given TR. As we increase our echo train length, we theoretically decrease our scan time and increase our SNR. The disadvantage of this parameter is that as we increase our echo train length, we increase the number of low-quality	Number

					echoes collected which lead to image blurring. Echo train length also promotes T2 contrast when increased.	
Gradient mode	GM	–	–	Use of grading	This parameter will change the way we ramp up our gradients. In other words, it changes our rise time rate. By doing this, we can potentially reduce artifacts in our image, reduced the loud noises produced from the MRI unit, and increase our echo spacing.	Drop Down Fast Normal
RF Mode	RFM	–	–	Use of grading	This parameter will change the way we perform our RF pulse. By increasing the duration of our RF pulse, we can decrease our SAR. We call this a soft pulse. This however can lead to increased scan times. On the contrary, if we decrease the duration of our RF pulse, we deliver a short burst of energy which can increase SAR.	Drop Down Fast Normal Low SAR
Normal Mode	NM	–	–	Use for SAR	This is a mode of operation that controls the time-varying magnetic fields (RF and gradient coils). It is a mode that is meant to not cause physiological stress to our patient.	Default
First Level	FL	–	–	Use for SAR	This mode of operation also controls the time-varying magnetic fields (RF and gradient coils). It is a mode that may allow for physiological stress to be presented to our patient. Medical supervision is required to go into this mode.	Radial button
Pixel shift	PS	[–	Use for grading	This describes how well we are able to represent chemical shift in our image. It defines how many pixels are needed to represent chemical shift. This is determined by field strength, receiving bandwidth, and frequency encoding.	Show Only
B-Value	BV	0	3000	Use for gradient	This parameter is unique to DWI. It determines how much diffusion we will see in our image. As we increase this parameter, we may will see more diffusion in our tissue. This will also reduce our SNR.	Drop Down 0 50 800 2000
VENC	VENC	0	500	Use for gradient	This parameter is unique to the phase contrast sequence. It determines the velocity of tissue we wish to be sensitive to. As we increase our VENC, we become more sensitive to faster moving tissue. It is important to identify the appropriate velocity of the tissue being imaged in order to accurately represent it.	Number field
Breath hold	BH	–	–	Use for calculation	When wishing to image our patient while they hold their breath, we need to use this parameter. It will allow us to image our patient over one or more breath holds. The number of acquisitions is a parameter that controls this.	Radial

High level parameters		
Parameter	Abbreviate	Calculate
SNR		$\frac{(((((FOV * 1.71) * (1/PE * 1.71) * (1/FE * 1.71) * (Slice\ thickness * 1.71) * Half\ Fourier * RectFOV * (1/RBW) * (NSA * 1.71) * (1 + OS)) / 2.5) / PI) * 1000 * Field\ Strength}}{((144 - (ETL + (Slice\ Thickness * 2.5) + (((1/RBW) * 1000) * 9)) - ((Half\ Fourier * 100 - 100) / 2))}$
CNR		$(144 - (ETL + (Slice\ Thickness * 2.5) + (((1/RBW) * 1000) * 9)) - ((Half\ Fourier * 100 - 100) / 2))$
Resolution		$(FOV/PE) \times (FOV/FE) \times Slice\ thick$
Scan time		$\frac{(((FS * (TR/ETL)) + TR) * (PE * averaged * half\ Fourier * rectfov) * (1 + phase\ oversampling))}{(ETL * Accelerate\ factor)}$
Temporal Resolution		TR*Segments Only Show when Retrospective button in the cardiac tab is selected
SAR		$\frac{Field\ strength^2 * B1+rms * (WARP) * RECTFOV + (oversampling))}{(ES) / 15}$ Popup a prompt: “SAR too High” If Normal Mode = (See SAR Table) First Level = (See SAR Table)
SED		Simulation 1: (SAR * Scan Time in sec.) + Simulation 2: (SAR * Scan Time in sec.) + Simulation 3: (SAR * Scan Time in sec.)

		+ Etc.
B1+rms		$\text{Flip angle}^2 * (((\text{Slices} * (\text{PE} * \text{Half Fourier})) * \text{RBW} / 100) / (\text{TR}))) / 100000$
PNS		

If SAR limit is met (See SAR Table), Show Pop up “SAR Limit Met” and give option to go to First Level Mode or adjust calculations again (retry)

SAR Table		
Coil	Normal Mode Limit	First Level Limit
Body RF Coil	2 W/kg	4 W/kg
Knee T-R	20 W/kg	40 W/kg
Brain R-only	3.2 W/kg	3.2 W/kg
Brain T-R	10 W/kg	20 W/kg
Surface R only	2 W/kg	4 W/kg
Foot R only	2 W/kg	4 W/kg
Breast R only	2 W/kg	4 W/kg
Shoulder R only	2 W/kg	4 W/kg

Template layout		
RITE Advantage		
Viewer	Viewer	Viewer
SAR: ____ W/kg SED: ____ J/kg B1+rms: ____ μT		
Scan Time: ____ Resolution: ____ x ____ x ____ SNR: ____		
<div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">Protocol</div> <div style="border: 1px solid black; padding: 2px;"> Simulation 1 Simulation 2 Simulation 3 Question 1 Simulation 4 <div style="text-align: right; font-size: small;">Total Scan Time: ____</div> </div> <div style="border: 1px solid black; height: 80px; margin-top: 10px; text-align: center; color: orange; font-weight: bold; font-size: 1.2em;">(Lesson Area)</div> <div style="background-color: green; color: white; text-align: center; padding: 5px; margin-top: 10px; width: 100px; float: left;">SCAN</div>	<div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">Contrast</div> <div style="border: 1px solid black; padding: 2px;"> TR: ____ ms Min. Max. TE: ____ ms Min. Max. Flip Angle: ____ degrees Saturation: <div style="border: 1px solid black; padding: 2px; display: inline-block;">Fat Water</div> Inverse Time: <div style="border: 1px solid black; padding: 2px; display: inline-block;">160 2500</div> Echo Train Length: ____ Regain Long. Mag. <input type="checkbox"/> </div> <div style="text-align: center; border: 1px solid black; padding: 2px; margin-top: 5px;">Advanced</div>	<div style="border: 1px solid black; padding: 2px; margin-bottom: 5px;">Resolution</div> <div style="border: 1px solid black; padding: 2px;"> # of Slices: ____ Slice Thickness: ____ mm Slice Gap: ____ mm FOV: ____ mm Image Matrix Frequency Encoding: <div style="border: 1px solid black; padding: 2px; display: inline-block;">512 448...</div> Phase Encoding: <div style="border: 1px solid black; padding: 2px; display: inline-block;">512 448...</div> </div> <div style="text-align: center; border: 1px solid black; padding: 2px; margin-top: 5px;">Misc</div>
Cardiac		

Viewer

Viewer

Viewer

SAR: ____ W/kg
SED: ____ J/kg
B1+rms: ____ μ T

Scan Time: ____ Resolution: __ x __ x __ SNR: ____

Protocol

Simulation 1
Simulation 2
Simulation 3
Question 1
Simulation 4

Total Scan Time: ____

(Lesson Area)

SCAN

Contrast

Advanced

Phase Encoding Direction:

A-P
I-R

Oversampling: ____ %

Number Averages: ____

Receiving Bandwidth: ____ Hz/px

Parallel Imaging: ☐

Half Fourier

100%
87% ..

RectFOV

100%
90%...

Confirm Frequency: ☐

Resolution

Misc

Echo Spacing: ____ ms

RF Mode:

Normal
Low SAR

Gradient Mode:

Fast
Normal

Normal Mode ☐ First Level ☐

MARS ☐

Measurement: ____

Temporal Resolution: ____

Breath holds: ☐

Cardiac

Viewer

Viewer

Viewer

SAR: ____ W/kg
SED: ____ J/kg
B1+rms: ____ μ T

Scan Time: ____ Resolution: __ x __ x __ SNR: ____

Protocol

Simulation 1
Simulation 2
Simulation 3
Question 1
Simulation 4

Total Scan Time: ____

(Lesson Area)

SCAN

Contrast

Advanced

Resolution

Misc

Retrospective
Trigger Window: ____
Trigger Delay: ____
Acquisition Window: +/- ____
Temporal Resolution: ____

Prospective
Segments: ____
Calc. Phase: ____


R-R

Trigger Delay

Acquisition Window

Trigger Delay

Cardiac



Part 3: Grading

Based on the data the user entered grades will be presented to them. The grade scale is:



Fail	Average	Pass
0-69%	70-85%	86-100%

This will be based on the following.


Parameter	Calculate	Fail Message	Average Message	Pass Message	Filter
SNR	= SNR Percent If pulse sequence all except spin echo and TR*2/FA = Less than 1, subtract 20% SNR	Your SNR was not sufficient enough to produce a diagnostic image.	Your SNR was acceptable	Great SNR	Noise filter = %
CNR	= CNR Percent	Your CNR was not sufficient enough to produce a diagnostic image.	Your CNR was acceptable	Great CNS	Gaussian filter = %
Resolution	= pixel size (FOV/PE x FOV/FE) PE = 100%< ± 3 default > 0% FE = 100%< ± 3 default > 0%	Your resolution was not sufficient enough to produce a diagnostic image.	Your resolution was acceptable	Great resolution	Sharpness filter = %
Breath hold	If breath hold button is selected and breath hold time is >30s=Fail 0% Breath hold time 20-29s=Average 80% Breath hold time <19s=pass 100%	Your breath hold was too long. Patient could not hold their breath for the entire sequence	Your breath hold was acceptable for the average patient.	Your breath hold was great and tolerable for most patients.	
Temporal Resolution	If retrospective is selected in the cardiac tab, show grade If >50ms Fail 0% If 40-49ms=Average 80% If <39ms=pass 100%	Your temporal resolution was too long. This will poorly represent moving tissues		Your temporal resolution was great	
Image Contrast	Uses Parameter vs Admin Selected	You have created the incorrect image contrast			= Assigned photo
Image Size	Filter only	Filter only			Zoom = User FOV / Admin FOV
Scan Time	< default input = pass	Your scan time was too long for this exam.	Your scan time was acceptable	Excellent scan time	
SAR	< default input when apply = pass	The patient exceeded their RF limit		SAR Limit Met	
SED	<limit when apply = pass	The patient exceeded their total RF dose			
B1+rms	< default input when apply = pass	The patient exceeded the B1+rms limit		B1+rms Limit Met	
Artifacts	= Average of Chemical shift, aliasing, partial volume				
<i>Chemical shift</i>	=Pixel shift If pixel shift is > 2.3, user fails. Otherwise they pass	Your technique presented a large pixel shift which lead to chemical shift artifact. Given your field strength, you should have used a larger receiving bandwidth or increased your frequency encoding in your image matrix			
<i>Aliasing</i>	Phase direction=RectFOV+FOV+Aliasing If phase direction (formula above) < box 2 size, use fails. Otherwise they pass	You have not taken consideration of tissue outside your field of view in the phase direction. You presented aliasing artifact			
<i>Partial Volume</i>	If slice thickness is > 8, user fails. Otherwise they pass	Your slice thickness was too thick and partial volume artifact was seen.			

Image Filtering


Stock Photo (image contrast user performed)




Noise Filter



Gaussian Filter



Sharp Filter



Zoom Filter

Image contrast will also be graded: If the selected parameters by the user = the image contrast defined in the builder, then this area is correct. The table below define the correct parameter for each given image contrast.

The image displayed and filtered to the user should be the image contrast they created.

If Spin echo					
Parameters	T1	T2	PD	IR	Hide Parameters
TR	400-900ms	>2000	>2000	STIR = >5000 FLAIR = >9000	VENC B value
TE	<30	>80	<30	-	
Flip Angle	>120	>120	>120	>120	
TI				1.5T TI = = STIR T1= = FLAIR 3.0T T1 = 220 = STIR T1 = 2500 = FLAIR	

If Gradient Echo				
Parameters	T1	T2	PD	Hide Parameters
TR	<250	>251	>251	ETL VENC B value
TE	<20	>21	<20	
Flip Angle	>45	<45	<45	

If DWI					
Parameters	EPI	Low Diffusion	Medium Diffusion	High Diffusion	Hide Parameters
TR	>3000	>3000	>3000	>3000	Change name ETL to Filling Factor VENC Flip Angle
TE	<200	<200	<200	<200	
B Value	0	<100	101-700	>701	

If Phase Contrast				
Parameters	Too Low	Just Right	Too High	Hide Parameters
TR	<10000	<10000	<10000	B Value Half Fourier RectFOV ETL
TE	<50	<50	<50	
Flip Angle	>10<90	>10<90	>10<90	
Venc	<set amount	Set amount	>set amount	

If SSFP		
Parameters	T2	Hide Parameters
TR	<20	B value Half Fourier VENC ETL
TE	<10	
Flip Angle	>20	

If TOF		
Parameters	Good	Hide Parameters
TR	<50	VENC B Value ETL
TE	<10	
Flip Angle	>5	

Grade Display

After every simulation submitted the user would get a popup of a score card.

- Grades should be set up to be exported
- Image result should be saved with images for (1 year).

Frontend

After the user clicks “scan” on the controls, a grade card of that sequence will be presented as a popup.

The image shows a frontend UI mockup. On the left is a control panel with three viewer windows at the top. Below them are input fields for SAR (W/kg), SED (J/kg), and B1+rms (μT), followed by Scan Time, Resolution (x, y, z), and SNR. A 'Protocol' dropdown menu lists Simulation 1, Simulation 2, Simulation 3, Question 1, and Simulation 4. Below the protocol is a 'Lesson Area' box. A 'SCAN' button is at the bottom left. An orange arrow points from the 'SCAN' button to a grade card popup on the right. The grade card popup has a 'Viewer' window at the top. Below it are input fields for SAR, SED, and B1+rms, followed by a 'Protocol' dropdown. The main area of the popup contains a 'Filtered Image' placeholder, a 'Total Score' field, and a 'Comments' text area. A 'SCAN' button is at the bottom of the popup.

When the activity is complete, a total score card will be demonstrated and saved under their dashboard. “Activity Results”

The image shows an 'Activity Results' UI mockup. At the top, it displays 'Activity Name' and 'Total Score'. Below this are three tabs: 'Sim 1', 'Sim 2', and 'Sim 3...'. An orange arrow points from the 'Sim 3...' tab to a 'Click to view grade card for specific simulation' text. Below the tabs is a list of parameters: SNR, CNR, Resolution, Scan Time, Temp Resolution, Image Contrast, SAR, SED, B1+rms, PNS, and Artifacts. A 'Comments' text area is below the list. A green 'Finish' button is at the bottom right.

Future Plans

1. Export Grades to Learndash
2. Dashboard for simulators: Show Completed activities and score
3. Free Play: have user control some of the conditions that build and exercise (Template, Field Strength, Pulse sequence, Image contrast, SAR limits, B1+rms limits, Scan Time limits)
4. Have a list of required activities that need to be completed with a pass score
5. Parameters displayed on different templates. (Template layout coming later) (Vendor to Vendor)
6. Multiple template view at once with breadcrumb to parameter.

Vendor to Vendor Example

Viewer

SAR: ___ W/kg
SED: ___ J/kg
B1+rms: ___ µT

Protocol

- Simulation 1
- Simulation 2
- Simulation 3
- Question 1
- Simulation 4

Total Scan time: ___

(Lesson Area)

SCAN

Scan Time: ___ Resolution: ___ x ___ x ___ SNR: ___

Contrast

Advanced

Phase Encoding Direction: A-P

Overampling: ___ %

Number Averages: ___

Receiving Bandwidth: ___ Hz/px

Parallel Imaging: ☐

Half Fourier: 100%
87% ..

RectFOV: 500%
47% ..

Confirm Frequency: ☐

Resolution

Misc

Echo Spacing: ___ ms

RF Mode: ☐ Normal ☐ Fast

Gradient Mode: ☐ Fast ☐ Normal

Normal Mode ☐ First Level ☐

MARS ☐

Measurement: ___

Temporal Resolution: ___

Breath holder: ☐

Cardiac

Contrast

Resolution

Geometry

Inline

Sequence

TR: ___

TE: ___

RBW: ___

ETL: ___

Template 1



Breadcrumb functionality

When a user selects a parameter in template 1, it will be found in template 2.

TO

Template 2



Terminology changes between templates (Send Filled Out Later)

Parameters	SIEMENS	GE	PHILIPS	Toshiba
TR				
TE				
TI				
Flip Angle				
Slice Thickness				
Slice Gap				
Field of View				
Phase Encoding				
Frequency Encoding				
Number of Slices				
Phase Direction				
Anti-Aliasing				

RectFOV				
Number of Signal Averages				
Number of Acquisitions				
Saturation				
Restore				
Half Fourier				
Parallel imaging				
Confirm Frequency				
Cardiac Gating				
Segments				
Calculated Phases				
Trigger window				
Trigger Delay				
Acquisition window				
Measurements				
Receiving Bandwidth				
Echo Spacing				
Echo Train Length				
Gradient mode				
RF Mode				
Normal Mode				
First Level				
Pixel shift				
B-Value				
VENC				
Breath hold				